


“Development application and single use probes for building Design Space of live cell concentration by dielectric spectroscopy”

BY: JOHN CARVELL

ABER INSTRUMENTS, Wales, UK



ABER Instruments Ltd – originally a spin out company from the University of **ABER**ystwyth, Wales.



Wales also has the longest station name in the world!



Fortunately it is shortened to Llanfair PG

Scope of Talk

- Principle of Dielectric Spectroscopy (DS) –emphasis on scanning
- Correlation with offline readings
- Case study Gedeon Richter
- Case study GSK
- Case study Fujifilm Diosynth
- Single use probe applications
- Conclusions



Industry Initiative

“Systems that promote greater product and process understanding can provide a high assurance of quality on every batch and provide alternative, effective mechanisms to demonstrate validation”

FDA Guidance for Industry: PAT—A Framework for Innovative Pharmaceutical Manufacturing and Quality Assurance, 2004.

Four types of tools for generation and application of process understanding:

1. **Multivariate tools for design, data acquisition and analysis**
2. **Process Analyzers**
3. **Process Control tools**
4. **Continuous improvements and knowledge management tools**


Guidance for Industry PAT — A Framework for Innovative Pharmaceutical Development, Manufacturing, and Quality Assurance

U.S. Department of Health and Human Services
Food and Drug Administration
Center for Drug Evaluation and Research (CDER)
Center for Veterinary Medicine (CVM)
Office of Regulatory Affairs (ORA)

Pharmaceutical CGMPs
September 2004



Principle of Dielectric Spectroscopy (DS)

- When an alternating electrical field is applied to an ionic solution of cells, a charge separation or polarization effect develops at the cell membrane
- Cells behave like a capacitor under an electric field
- The magnitude of this charge separation can then be measured by the capacitance of the solution at varying radio-frequencies in the range 0.1-20MHz
- The resulting capacitance is generally proportional to the total enclosed bio-volume, which itself is directly correlated with cell density or cell weight



Dielectric Spectroscopy (DS) in bioreactors

- Dielectric Spectroscopy Key Outputs:**
- Capacitance & Dielectric Data Over Time
 - Trends about Cell Viability Density




Graphics courtesy of Bend Research, Oregon



Correlation of online capacitance data with offline readings



- Plot of capacitance at specific single frequency versus VCD should produce straight line correlation with an R^2 value close to 1
- Correlation during log phase of cell growth is well maintained
- Divergence observed as cells enter stationary/death phase at which point cell viability decreases and LDH levels increase proportionally
- Published data on models for correcting this divergence available *eg* Bend, UMASS, University of Manitoba




CHO Cell Density Monitoring in Bioreactor Culture – Comparison of four methods (data courtesy) University of Manitoba)




Fed-Batch HEK cells (data courtesy of CNRC, Canada)




0.6 MHz



2.2 MHz




10 MHz



Change in shape of Spectrum- fundamentals of frequency scanning


Control

Timeseries of Normalized Capacitance vs. Frequency



CHO cell undergoing induced apoptosis

Timeseries of Normalized Capacitance vs. Frequency



Application of dielectric spectroscopy for monitoring high cell density in monoclonal antibody producing CHO cell cultivations

László Párta · Dénes Zalai · Sándor Borbély ·
Ákos Putics

Case Study 1

Gedeon Richter, Hungary



GEDEON RICHTER

Contact Dénes Zalai
zalaid@richter.hu



Details of study at Gedeon Richter


- 6 runs with 1L bioreactors and ATF2 technology (Refine,USA) with CHO expressing Mab
 - High cell densities up to 75×10^6 cells/ml
 - Utilized Aber Biomass Monitor 220 with AberScan software
 - Scans at 25 pre-defined frequencies every 8 minutes
 - Spectral data exported for linear regression or multivariate data analysis
-
- Cells/ml derived from Aber using 3 different methods:
 - Simple measurement at 580KHz (linear modelling)
 - Cole- Cole
 - PLS using Umetrics SIMCA software



GEDEON RICHTER



Linear Modelling at 580KHz



- at-line measured VCD values throughout Run 1-4 were plotted against predicted VCD values and depicted as triangles
- linear equation for prediction:
$$\text{VCD} = 414492 \times C_{580} + 1.187 \times 10^6$$
- linear model over-estimated VCD at the later phases of the cultivation




GEDEON RICHTER

ABER
TRUSTED TECHNOLOGY

Deriving VCD using Cole-Cole using frequency scanning between 100 and 20,000KHz

Theoretical equations for β -dispersion curve (Cole-Cole equation)




$$\Delta\epsilon = \frac{9pRC_m}{4(1+RG_m(1/\sigma_i+1/2\sigma_a))^2} \approx \frac{9pRC_m}{4} = 3\pi R^4 C_m \cdot VCD$$

$$f_c = \frac{2\sigma_i\sigma_a + RG_m(\sigma_i + \sigma_a)}{2\pi RC_m(\sigma_i + 2\sigma_a)} \approx \frac{1}{2\pi RC_m(1/\sigma_i + 1/2\sigma_a)}$$

- ❖ The known values of C_m , σ_i , and σ_a together with $\Delta\epsilon$ and $f_c \rightarrow VCD$ and R (cell radius) can be calculated from the above two equations
- ❖ The known values of C_m , σ_i , and σ_a together with $\Delta\epsilon$ and $f_c \rightarrow$ Assume that cell radius, R is constant during culture \rightarrow

$$\Delta\epsilon \propto pR \propto VCD \cdot R^4 \propto VCD$$


Off-line cell counts (VCD) for 6 runs with REFINE system



GEDEON RICHTER


ABER
TRUSTED TECHNOLOGY

Score plots for 6 runs




- Multivariate analysis (PCA) of capacitance spectra resulted in detection of measurement disturbances
- Principle component $t(1)$ captured 92% of variance in spectra and showed same trend for all 6 runs
- **DIELECTRIC SPECTRA NOT INFLUENCED BY BATCH TO BATCH DIFFERENCES OF CRITICAL CULTIVATION PARAMETERS INDICATING ROBUSTNESS OF probe**
- VCD prediction with multivariate PLS model developed

CHO with continuous feed-Comparison of 3 methods to predict VCD



CHO with continuous feed-Comparison of 3 methods to predict VCD



Case study 2

GSK, Stevenage, UK

Dielectric Spectroscopy and its application in process development

Andrew Heinrich

Biopharm R&D, GSK Stevenage


PAT & QbD Forum, Feb 18th

Goettingen, Germany





Use of capacitance in scale-up and technology transfer

- Background
 - Cytokine produced in modified CHO cells
 - Initial cell line selection carried out in shake flask studies
 - Process scoping in 2L bioreactor system
- In-situ Aber FUTURA Probe used in 2L glass bioreactors with multiple frequency scanning
- What information can we glean from this data?
 - Relatively similar process profiles
 - Perturbations in capacitance profiles appear to correspond with feed points in process
 - Can this data be used to help build a 'process model' for scale up (and scale-down)?



- GSK developed an XL-fit macro for rapid and reproducible automatic analysis of scanning data sets


			XLFIT	E _{max}	Slope	EC50	EO	Delta E _{max}	Biovolume P
12139	15582	20000							
-0.05384	-0.02003	-0.02867019	#Ok	-0.02213	67.9377	50	0	-0.02212655	-1.10632767
-0.05609	-0.02926	-0.02333686	#Ok	0.15623	-84.3199	50	0	0.1562311	7.81155525
-0.06892	-0.04473	-0.04530769	#Ok	0.12267	-10.0714	1096.66	0	0.12266904	134.526142
-0.05755	-0.0459	-0.03802019	#Ok	0.24021	-6.3871	1017.35	0	0.24021494	244.383568
-0.06545	-0.03558	-0.03590353	#Ok	0.30406	-6.23227	1057.29	0	0.30405766	321.476758
-0.07727	-0.03394	-0.03028269	#Ok	0.40628	-5.25257	1053.34	0	0.40628334	427.955626
-0.07307	-0.04495	-0.02983686	#Ok	0.52413	-4.33303	973.588	0	0.52412527	510.281922
-0.0644	-0.04436	-0.03675353	#Ok	0.65062	-4.10737	954.47	0	0.65062089	620.99819
-0.06359	-0.04909	-0.03325353	#Ok	0.85102	-3.63535	910.199	0	0.85102195	774.599067
-0.03963	-0.01646	0.005313141	#Ok	1.09639	-3.23762	869.539	0	1.09638907	953.352867
-0.04287	0.00657	0.029250641	#Ok	1.30887	-3.082	848.211	0	1.30886902	1110.1977
-0.04326	-0.00259	0.018208974	#Ok	1.52516	-2.88532	825.095	0	1.52515804	1258.40068
-0.04691	-0.01544	0.013429808	#Ok	1.73297	-2.77729	813.643	0	1.73296655	1410.01534
-0.0404	-0.01336	0.022763141	#Ok	1.99487	-2.64025	781.499	0	1.99486599	1558.98574
-0.0071	0.02674	0.055717204	#Ok	2.31913	-2.41333	760.452	0	2.31912931	1763.58753
0.01036	0.05772	0.096129808	#Ok	2.60404	-2.33824	759.638	0.01012	2.59391923	1970.44014
0.00974	0.05008	0.087875641	#Ok	2.81761	-2.33608	762.594	0.01426	2.80334183	2137.81106
0.02664	0.06881	0.095208974	#Ok	3.07825	-2.28322	756.015	0.01675	3.06149186	2314.53463
0.02434	0.07905	0.105646474	#Ok	3.40322	-2.16663	730.877	0.02845	3.37476971	2466.54145
0.01963	0.06479	0.091971474	#Ok	3.57071	-2.21813	747.204	0.02232	3.54839176	2651.37382
0.00873	0.05317	0.080313141	#Ok	3.8613	-2.16156	738.955	0.01006	3.85123908	2845.89176
-0.023	0.00977	0.026117308	#Ok	3.90199	-2.20813	731.644	0	3.90198695	2854.8648
-0.03276	0.00119	0.016242308	#Ok	4.27703	-2.12814	707.708	0	4.27703031	3026.88694
-0.07777	-0.06187	-0.05085769	#Ok	4.35397	-2.24076	712.702	0	4.3539736	3103.08611



Use of capacitance in scale-up and tech transfer

- Analysis of various scanning data parameters highlights clearer differences between some of the processes
- Using this information one can eliminate runs that appeared previously to be consistent and therefore make more informed conclusions of the 'predictability' of the process conditions
- Further processing of data is still required to fully analyse these results and gain a greater confidence in the 'process model'


Biovolume profiles for Bioreactor runs



Multivariate analysis



- PCA modelling can then be used to further interpret the key components and specific characteristics of individual batches
- Pre-processing of the raw data is critical in this sort of analysis in order to ensure that inferences from the statistical tools are 'genuine'
- This analysis then offers a method for investigating design space and process 'limits' suitable for development and scale up



Using DS in scale up/scale down experiments

Using conventional probes volumes down to 100ml can be monitored




Optimization of human pluripotent stem cell suspension cultures in controlled, stirred bioreactors

Ruth Olmer¹, Ulrich Martin¹ Robert Zweigerdt¹

*Leibniz Research Laboratories for Biotechnology and Artificial Organs (LEBAO),
REBIRTH Cluster of Excellence, Hannover Medical School, Carl-Neuberg-Str. 1, 30625 Hannover, Germany*


Samples down to 1ml can be taken from shake flasks for DS using an inverted probe



Scale up data from GSK

- Capacitance as a scale parameter from shake flask to 50L single use bioreactor
- Inverted probe technology from Aber Instruments provides real-time offline capacitance readings that correlate well with online data


Correlation of Off-line v *In-situ* Capacitance probes at scale




▶ Aims : Acquire capacitance spectra data & compare conventional calibration to PLS calibration

- Glass 3L bioreactor
- CHO DG44 clone
- Fed batch process
 - Cell count & Viability data from Vicell
 - Conventional calibration carried out after run and applied to historical capacitance data
 - PLS calibration carried out after run and applied to historical capacitance data

Case study 3
Fujifilm Diosynth, UK



Conventional Calibration showing online cell diameter from ABER SCADA software




FUJIFILM

Diosynth
biotechnologies

ABER
TRUSTED TECHNOLOGY

Principal Component Analysis


Raw scores plot of B-Dispersion Spectra



- PCA analysis of Raw B-dispersion Spectra from CHO Mammalian cell culture
- P1 – Strong Biomass concentration signal
- PC2 and 3 Remaining Orthogonal variance in signal

VCD- Variable importance plots vs Frequency

131216 PLS predictionset.M2 (PLS), Viable cells only
 VIP[Comp. 3]




For VCD the lower frequencies correlate more strongly than the higher frequencies as expected

Viability-Variable importance plots vs Frequency

-Two peaks in the VIP plot vs frequency with a peak in the mid range.

-change in the shape of the peak in the mid range over time probably due to a shift in the cell size distribution.


131216 PLS predictionset.M3 (PLS), % Viability Only
 VIP[Comp. 6]



PLS Calibration

FUJIFILM

Diosynth
biotechnologies




ABER
TRUSTED TECHNOLOGY

Single use applications





Results from Aber single use probe in a Sartorius rocking bag



Flow past cells for use on external loops with silicone tubing



First version of Flow past cell will be tested at Shire this month



Conclusions

- Spectral scanning of capacitance becoming increasingly used in PAT
- Software based on Cole-Cole or PLS can be used to correct VCD during death phase ,derive viability and in scale up studies
- Promising technology for studies on apoptosis on certain cell lines
- Development of single use probes has enabled technology to be applied to single use bioreactors

