



Urological

pAguamedicina Structural Hydrogel
Medical Device Material

Scott M. Epstein

CEO / CTO, Q Urological & Boston Scott Corporations

A wholly owned subsidiary of Boston Scott Corporation

polyaguamedicina used as a **medical device** (Catheter) raw material should provide an immediate reduction in health care costs; increased patient comfort and quality of care.



polyaguamedicina is a hydrolyzed polyacrylonitrile; used in a unique proprietary manufacturing process produces structural 100% hydrogel devices, comprised mostly of water (or aqueous media) essentially built out of water.

polyaguamedicina is currently **FDA approved** (for pediatric ureteral stents), exhibiting true inherent **bio-compatibility** and mechanical **compliance** is expected to reduce the incidence of **encrustation; infection & tissue irritation**.

polyaguamedicina provides antifouling rather than antimicrobial performance; that is no additives are required whereby biofilms are not expected to adhere and migrate

Comparison of Medical Device (Catheter) Materials

This is not a hydrogel coating, the hydrogel material is the catheter; there is no “coating” to erode.



Conventional thermoplastic catheter material example:
100% polymer



Conventional thermoplastic catheter material example:
80% polymer / 20% radio-pacifier



Hydrophilic polyurethane example:
60% polymer, 20% radio-pacifier,
20% water



pAguamedicina example:
10% polymer, 90% water
Radio-pacifier, negligible



Typical urological short term implant complication




Polyurethane sample 10 X magnification of surface
(typical encrustation)

Bio-compatibility (In-vitro) of catheter materials

Potential to eliminate bio-film formation & encrustation.

7 days submersion in human urine at 98°F

Typically competitive sample



pAguamedicina stent sample




10 X magnification of
surface (no sign of
encrustation)

Bio-compatibility (continued)


<u>ISO / FDA / CE requirements</u>		
<u>Test description / name</u>	<u>ISO ref. #</u>	<u>Score</u>
Cytotoxicity using 1X Minimal Essential media.	10993: 5	<u>PASS:</u> No evidence of cell lysis or toxicity was observed.
Muscle implantation (4 week) study with Histo-pathology	10993: 6	<u>PASS:</u>
Intracutaneous Study, 0.9% Sodium Chloride Solution Extract	10993: 10	<u>PASS:</u> There was no erythema or no edema from extract tested injected intracutaneously into rabbits.
Intracutaneous Study, Sesame Oil NF Extract		
Systemic Toxicity Study; 0.9% Sodium Chloride Solution Extract	10993: 11	<u>PASS:</u> There was no mortality or evidence of systemic toxicity with either extract.
Systemic Toxicity Study; Sesame Oil NF Extract		
Genotoxicity, Bacterial Reverse Mutation Study, 0.9% Sodium Chloride Solution Extract	10993: 3	<u>PASS:</u> Considered to be non-mutagenic
Genotoxicity, Mouse Lymphoma Assay, 0.9% Sodium Chloride Solution Extract		<u>PASS:</u>
<u>Analytical testing</u>		
Exhaustive Extraction , Lipophilic; Ethanol Extract; Non volatile residue		<u>PASS</u> <u>All results were acceptable. (PPB levels of detectability)</u>
Inductively Coupled Plasma Spectroscopy , Full Scan- Extractable Metals		
Determination of Extractable Semi-Volatile Organic Compound by GC / MS		
LC/MS Chromatographic Screen for Extracts		

Radiopacity

(Can also be engineered for visibility in CT, MR and Ultrasound environments)




Mechanical Characteristics



Comparative column strength


Push-ability




Anchorage details

Conventional pigtail anchorage end;

whereby loop diameter * 3.14 = anchorage compliance; essentially the unwinding or the loop due to peristaltic and or anatomic movement. However testing indicates that at about 50% of circumferential travel; pigtail becomes hung-up, or stuck, contributing we believe to pain and or patient **discomfort**.



Drainage, must conduct through the lumen entering in through ports and or the open end of the stent in the kidney. However especially small tubes (< 4F) ports are ineffective and easily clogged.



As shown, assumes patient is vertical, when resting flow into stent would be along its center axis.

Although they can be provided, pAguamedicina stents will not require ports; rather flow will be conducted form the wide open bulbous end. Similarly due to the very hydrophilic nature of the stent material; urine will be absorbed and transported accordingly into through the stent

Dynamic analysis




Package

Product Design Opportunity

Shapes, profiles, geometry, Hybrids

(Can incorporate shape memory, be used as a unique device material or combined as a hybrid, with other materials or components)

Currently approved: 510K cleared 01/20/2010 **pAguamedicina™** Structural Hydrogel Pediatric Ureteral Stent,



Summary of Expectations

Very biocompatible medical device material

Instantly slippery when saturated with saline, water or bodily fluids

Reduced coefficients of friction when sliding coaxially over or against dissimilar materials.

Can be saturated with almost any aqueous media, exhibiting delivery and absorption gradients; remaining stable and impervious to most solvents.

Saturated in aqueous solutions (pH: 5.5 to 7.0) exhibits (about 20A Shore durometer); however can be engineered to exhibit specific characteristics.

Can not ETO or steam sterilize; however Radiation sterilization method **does not** effect polymer characteristics.